

**E008: Auto-disable syringe for fixed dose immunization**

PQS code: E008/052
Description: AD syringe 0.05ml
Manufacturer's reference: OnejectII 0.05ml AD syringe
Manufactured in: Indonesia
Company: PT Oneject Indonesia
Address: Sentul Industrial Estate
Jl.Olympic Raya Kav B9
Bogor 16810,
Indonesia

Telephone: +62 21 879 17 422

Email: marketing@oneject.co.id

Web address: www.oneject.co.id

Specifications

Graduations:	0.05ml	Quality standard:	ISO 13485,
Syringe material(s):	Polypropylene	Markings:	CE mark
Fixed needle size:	27G x 3/8" (0.40 x 10 mm)	Pieces per carton:	1600
AD mechanism:	Locked / broken plunger	Volume per carton (m3):	0.051m3
AD location:	Upon completion of injection	Weight per carton (kg):	8.23kg
Number of components:	3 pieces	Minimum order:	16,000
Other needle options:	Yes	Incoterms:	FOB
Primary packaging:	Blister pack	Prices per unit:	<i>Upon request to manufacturer</i>
Year base price:	2013		

Comments:

Current PQS status: *pre-qualified ; 05 Jun 2013*

Valid until : *May 2021*

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased