




PERFORMANCE QUALITY SAFETY



E013: Hypodermic syringes with reuse prevention feature

PQS code: **E013/108**
Description: RUP Syringe 2 mL
Manufacturer's reference: Oneject II RUP Syringe 2 mL
Manufactured in: Indonesia
Company: PT Oneject Indonesia
Sentul Industrial Estate
Jl.Olympic Raya Kav B9
Bogor 16810,
Indonesia
Address:

Telephone: +62 21 879 17 422
Email: marketing@oneject.co.id
Web address: www.oneject.co.id

Specifications

Graduations:	0.1 mL	Needle size:	23G x 1 ¼"
Syringe material(s):	Polypropylene	Other needle options:	Yes
Number of components:	3 pieces	Needle fixation:	Fixed
RUP mechanism:	Locked / broken plunger	Primary packaging:	Paper blister pack
RUP location:	Upon completion of injection	CE mark:	Yes
RUP type:	Type I B		
Pieces per carton:	1600	Volume per carton (m3):	0.09 m3
Minimum order:	1600	Weight per carton (kg):	10.6 Kg
Year base price:	2015	Incoterms:	FOB
Price per unit:	Upon request to manufacturer		
Quality standard:	- ----- - ISO 9001 - ISO 13485 -		

Comments:

Current PQS status: *pre-qualified; 21 Sep 2015*

Valid until : *May, 2021*

Note: If Current PQS status is 'Suspended' or 'Withdrawn',