

**E013: Hypodermic syringes with reuse prevention feature**

PQS code: E013/019
Description: RUP syringe 3ml
Manufacturer's reference: Oneject K1 3ml
Manufactured in: Indonesia
Company: PT Oneject Indonesia
Sentul Industrial Estate
Jl.Olympic Raya Kav B9
Address: Bogor 16810,
Indonesia

Telephone: +62 21 879 17 422
Email: marketing@oneject.co.id
Web address: www.oneject.co.id

Specifications

Graduations:	0.1ml	Needle size:	23G x 1 1/4"
Syringe material(s):	Polypropylene	Other needle options:	Yes
Number of components:	3 pieces	Needle fixation:	Fixed
RUP mechanism:	Locked & broken plunger	Primary packaging:	Paper blister pack
RUP location:	Upon completion of injection	CE mark:	Yes
RUP type:	Type I B		
Pieces per carton:	1800	Volume per carton (m3):	0.11 m3
Minimum order:	18000	Weight per carton (kg):	13.6 kg
Year base price:	2010	Incoterms:	FOB
Price per unit:	On request to manufacturer		

Quality standard: - ISO 9001 - ISO 13485 -

Comments:

Current PQS status: *pre-qualified; 07 May 2007*

Valid until : *May, 2021*

Note: If Current PQS status is 'Suspended' or 'Withdrawn',