

**E008: Auto-disable syringe for fixed dose immunization**

PQS code: E008/030
Description: AD syringe 0.05 ml
Manufacturer's reference: K1 Oneject 0.05 ml
Manufactured in: Indonesia
Company: PT Oneject Indonesia
Address: Sentul Industrial Estate
Jl.Olympic Raya Kav B9
Bogor 16810,
Indonesia

Telephone: +62 21 879 17 422

Email: marketing@oneject.co.id

Web address: www.oneject.co.id

Specifications

Graduations:	0.05ml	Quality standard:	ISO 13485,
Syringe material(s):	Polypropylene	Markings:	CE 0044
Fixed needle size:	26G x 3/8" (0.45 x 10 mm)	Pieces per carton:	1600
AD mechanism:	Locked / broken plunger	Volume per carton (m3):	0.06 m3
AD location:	Upon completion of injection	Weight per carton (kg):	6.3 kg
Number of components:	3 pieces	Minimum order:	16,000
Other needle options:	Yes	Incoterms:	FOB
Primary packaging:	Blister pack	Prices per unit:	<i>On request</i>
Year base price:	2007		

Comments:

Current PQS status: *pre-qualified ; 21 May 2007*

Valid until : *May 2025*

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased