

**E013: Hypodermic syringes with reuse prevention feature**

**PQS code:** E013/019  
**Description:** RUP syringe 3ml  
**Manufacturer's reference:** Oneject K1 3ml  
**Manufactured in:** Indonesia  
**Company:** PT Oneject Indonesia  
Sentul Industrial Estate  
Jl.Olympic Raya Kav B9  
Bogor 16810,  
Indonesia

**Telephone:** +62 21 879 17 422  
**Email:** [marketing@oneject.co.id](mailto:marketing@oneject.co.id)  
**Web address:** [www.oneject.co.id](http://www.oneject.co.id)

**Specifications**

<b>Graduations:</b>	0.1ml	<b>Needle size:</b>	23G x 1 1/4"
<b>Syringe material(s):</b>	Polypropylene	<b>Other needle options:</b>	Yes
<b>Number of components:</b>	3 pieces	<b>Needle fixation:</b>	Fixed
<b>RUP mechanism:</b>	Locked & broken plunger	<b>Primary packaging:</b>	Paper blister pack
<b>RUP location:</b>	Upon completion of injection	<b>CE mark:</b>	Yes
<b>RUP type:</b>	Type I B		
<b>Pieces per carton:</b>	1800	<b>Volume per carton (m3):</b>	0.11 m3
<b>Minimum order:</b>	18000	<b>Weight per carton (kg):</b>	13.6 kg
<b>Year base price:</b>	2010	<b>Incoterms:</b>	FOB
<b>Price per unit:</b>	On request to manufacturer		

**Quality standard:** - ISO 9001 - ISO 13485 -

**Comments:**

**Current PQS status:** *pre-qualified; 07 May 2007*

**Valid until :** *May 2025*

Note: If Current PQS status is 'Suspended' or 'Withdrawn',