

**E013: Hypodermic syringes with reuse prevention feature**

PQS code: E013/024
Description: RUP syringe 5ml
Manufacturer's reference: Oneject K1 5ml
Manufactured in: Indonesia
Company: PT Oneject Indonesia
Sentul Industrial Estate
Jl.Olympic Raya Kav B9
Bogor 16810,
Indonesia

Telephone: +62 21 879 17 422
Email: marketing@oneject.co.id
Web address: www.oneject.co.id

Specifications

| | | | |
|------------------------------|------------------------------|--------------------------------|--------------------|
| Graduations: | 0.2ml | Needle size: | 21G x 1 1/2" |
| Syringe material(s): | Polypropylene | Other needle options: | Yes |
| Number of components: | 3 pieces | Needle fixation: | Detachable |
| RUP mechanism: | Plunger locked and broken | Primary packaging: | Paper blister pack |
| RUP location: | Upon completion of injection | CE mark: | Yes |
| RUP type: | Type I B | | |
| Pieces per carton: | 1200 | Volume per carton (m3): | 0.09m3 |
| Minimum order: | 12000 | Weight per carton (kg): | 11.1 kg |
| Year base price: | 2008 | Incoterms: | FOB |
| Price per unit: | On request to manufacturer | | |

Quality standard: - ISO 9001 - ISO 13485 -

Comments:

Current PQS status: *pre-qualified; 21 May 2007*

Valid until : *May 2025*

Note: If Current PQS status is 'Suspended' or 'Withdrawn',