

**E013: Hypodermic syringes with reuse prevention feature**

<b>PQS code:</b>	<b>E013/108</b>
<b>Description:</b>	RUP Syringe 2 mL
<b>Manufacturer's reference:</b>	Oneject II RUP Syringe 2 mL
<b>Manufactured in:</b>	Indonesia
<b>Company:</b>	PT Oneject Indonesia Sentul Industrial Estate Jl.Olympic Raya Kav B9 Bogor 16810, Indonesia
<b>Address:</b>	
<b>Telephone:</b>	+62 21 879 17 422
<b>Email:</b>	<a href="mailto:marketing@oneject.co.id">marketing@oneject.co.id</a>
<b>Web address:</b>	<a href="http://www.oneject.co.id">www.oneject.co.id</a>

**Specifications**

<b>Graduations:</b>	0.1 mL	<b>Needle size:</b>	23G x 1 ¼"
<b>Syringe material(s):</b>	Polypropylene	<b>Other needle options:</b>	Yes
<b>Number of components:</b>	3 pieces	<b>Needle fixation:</b>	Fixed
<b>RUP mechanism:</b>	Locked / broken plunger	<b>Primary packaging:</b>	Paper blister pack
<b>RUP location:</b>	Upon completion of injection	<b>CE mark:</b>	Yes
<b>RUP type:</b>	Type I B		
<b>Pieces per carton:</b>	1600	<b>Volume per carton (m3):</b>	0.09 m3
<b>Minimum order:</b>	1600	<b>Weight per carton (kg):</b>	10.6 Kg
<b>Year base price:</b>	2015	<b>Incoterms:</b>	FOB
<b>Price per unit:</b>	Upon request to manufacturer		
<b>Quality standard:</b>	- ----- - ISO 9001 - ISO 13485 -		

**Comments:**

**Current PQS status:** *pre-qualified; 21 Sep 2015* Valid until : *May 2025*

Note: If Current PQS status is 'Suspended' or 'Withdrawn',